

## Media Recipes

### **18% MGM (Modified Growth Media)**

Developed by Mike Dyall-Smith for growing *Haloferax volcanii* DS2

144 g NaCl  
18 g  $\text{MgCl}_2 \cdot 6\text{H}_2\text{O}$   
21 g  $\text{MgSO}_4 \cdot 7\text{H}_2\text{O}$   
4.2 g KCl  
8 ml 1 M Tris.Cl pH 7.5  
5 g Oxoid\* peptone  
1 g Yeast Extract\*  
to 1 L dH<sub>2</sub>O

You may have to gently heat the water to have the salts dissolve completely. Add salts slowly to a beaker with about 600 ml water and a stirbar. When everything is dissolved, bring volume up to 1 L.

If plates are desired, add 2% agar before autoclaving.

Autoclave 20 minutes. After autoclaving, add 6 ml of sterile 0.5 M  $\text{CaCl}_2 \cdot 2\text{H}_2\text{O}$

\*The Oxoid brand of peptone is required as the field has discovered inhibitory contaminants (bile salts) in brands such as Difco.

\*\*Recently, the Difco brand of yeast extract has been shown to inhibit *H. volcanii* growth, though the identity of the inhibitor is unknown. However, if you have Difco yeast extract from pre-2017, it will be fine.

### **CM (Complete Media)**

Hackett and Dassarma, 1989, for growing *Halobacterium salinarum* NRC-1

250 g NaCl  
20 g  $\text{MgSO}_4 \cdot 7\text{H}_2\text{O}$   
10 g Oxoid peptone  
3 g (Tri)sodium citrate  
2 g KCl  
to 1 L dH<sub>2</sub>O

Again, you may have to gently heat the water to have the salts dissolve completely. Add salts slowly to a beaker with about 600 ml water and a stirbar. When everything is dissolved, pH to 6.8 with NaOH and bring volume up to 1 L.

Add 2% agar for plates.

Autoclave 20 minutes.

We have many more media recipes in the Schmid lab if you have specific needs or trouble finding these components. Contact Cindy at [cynthia.darnell@duke.edu](mailto:cynthia.darnell@duke.edu) or [darnell.cynthia@gmail.com](mailto:darnell.cynthia@gmail.com) for discussion or more recipes.